



MOTIFMESHTM


High Performance Monolayer PTFE Mesh



Tailored for Intraperitoneal (IPOM) Repair



- ULTRA THIN FOR LAPAROSCOPIC DELIVERY
150 micron/ 0.006"

- NO ORIENTATION REQUIRED AS EQUAL STRENGTH IN EACH DIRECTION



- TRANSPARENT MONOLAYER MESH⁷ CAN BE PLACED ON EITHER SIDE

- SMOOTH BIOCOMPATIBLE PTFE MATERIAL
- MACROPOROUS FOR OPTIMAL TISSUE INTEGRATION
- INHERENT ANTI-ADHESIVE SURFACE PROPERTIES




Proven Performance Across Key Criteria

Performance Criteria	v. Flatsheet PTFE Mesh	v. Composite Mesh	v. Biologic Mesh
1: Mesh Infection Resistance	Higher	Equivalent - Lower	Equivalent - Lower
2: Tissue Incorporation	Higher	Equivalent - Higher	Higher
3: Visceral Adhesion	Lower	Equivalent - Lower	Higher
4: Shrinkage	Lower	Equivalent	Lower
5: Relative Pricepoint	ca. 60%	ca. 50% - 90%	ca. 10% - 30%

Product comparisons carried out in pre-clinical and clinical trials - Ref: 1,2,3,5,6


1. MESH INFECTION RESISTANCE



Reference 3 - Sanders et al, 2013

- MotifMESH condensed PTFE material structure inhibits bacterial formation better than competing products.^{1, 2, 3}

2. TISSUE INCORPORATION




Reference 2 - Deerenberg et al, 2012

- Study shows MotifMESH has faster Abdominal incorporation after 28 days than any other synthetic or biologic mesh.²
- Macroporous design promotes rapid healing and dense collagen formation.^{5, 6}

3. VISCERAL ADHESION


- Ultra smooth, non porous PTFE material provides inherent anti-adhesive qualities.⁶
- Image below showing good peritoneal covering with vascularisation of new peritoneum - homogeneous covering over undulating mesh with no adhesions to the mesh - one to a tack.

In-vivo study - MotifMESH 1yr post implantation



4. SHRINKAGE

- Shrinkage 10-15%, in line with most composite and non-composite mesh after 90 days, versus 23% shrinkage for biologic mesh.²
- Facilitates repair of large hernias.



References

1. "Macroporosity and hydropobicity of surgical meshes reduce in vivo staphylococcus aureus infection and anchorage" Voskerician et al. (Presented at 4th International Hernia Congress, Sept. 2009)
2. "Experimental Study on Synthetic and Biological Mesh Implantation in a Contaminated Environment", Deerenberg et al. (BJS Journal Vol. 99, Issue 12, Dec. 2012)
3. "An in vitro study assessing the effect of mesh morphology and suture fixation on bacterial adherence", D Sanders et al. (Springer Link, Hernia, Vol. 17, Issue 6, 779-789, Dec. 2013)
4. "Foreign body reaction to meshes used for the repair of abdominal wall hernias," Klinge, U. et al, Eur J Surg (1999); 165: 665-673
5. "Macroporous condensed poly(tetrafluoroethylene) I. In vivo inflammatory response and healing characteristics", Voskerician, G et al, J Bio Mat res (2006); 76A: 232-242.
6. "Macroporous condensed poly(tetrafluoroethylene). II. In vivo effect on adhesion formation and tissue integration", Voskerician, G et al, J Biomed Mater Res A. 2007 Aug; 82(2):426-35.)
7. In house benchmark study - data on file.